



**THE BILL**  
YOUR REGULAR DOSE OF EDUCATION  
FOR INFECTION PREVENTION

**FLUIDSHIELD\* SURGICAL N95 DUCKBILL RESPIRATOR  
FITS THE BILL  
QUALITY, COMFORT & PROTECTION**

**The HALYARD  
Education Podcast**

[LISTEN HERE](#)



**NOT ALL N95 RESPIRATORS OFFER THE SAME LEVEL OF PROTECTION**

With so many new suppliers in the global market and so much counterfeit product driven by demand, the recommendation is for healthcare workers not to assume that ANY mask does the job properly. When considering which N95 Respirator provides the right level of protection, healthcare workers should request from their mask supplier the testing and other necessary information.

## What questions should I ask my N95 Respirator supplier?

### Does the mask meet the standards?

Ensure your N95 Respirator meets one of the face mask standards either; AS4381:2015 (Australian), ASTM F2100-11 (North American), or EN14683 (European). These standards cover two key performance criteria:

- **Filtration:** how well the mask filters smaller microorganisms
- **Breathability:** how easy it is to breathe through the mask

### How effective is the mask against fluid/splashes?

If you are going to use the mask where the risk for fluid is high, ensure that you request information on how well the mask resists splashes. If you require a high level of splash resistance, ensure your mask is a Surgical N95 Respirator offering the highest splash resistance and rated to ASTM F2100-11 / AS4381:2015 Level 3.



### Is it NIOSH approved?

NIOSH is the National Institute for Occupational Safety and Health, a division of the U.S. Centers for Disease Control and Prevention. NIOSH is an independent testing body, globally recognised as the N95 specialist endorsing and validating N95 product performance versus respiratory standards. NIOSH provides assurance of product quality and tests respirators' ability to protect the wearer from dust and liquid droplets in the air.

Selecting an N95 Respirator with **NIOSH approval will provide the assurance that the mask will offer you the respiratory protection needed.**

### Is it approved by the TGA?

The Therapeutic Goods Administration (TGA) regulates N95 Respirators to ensure they are of a high quality, providing the protection they claim to. In order to be registered with the TGA, they must be entered in the ARTG. Request the ARTG number from your mask supplier to ensure **your mask has been approved by the TGA.**

**DON'T WEAR THE RISK.  
WEAR HALYARD.**

## HALYARD FLUIDSHIELD\* Surgical Level 3 N95 Respirators:

- ✓ **Meet the Standards** for single use face masks AS4381, ASTM F2100-11 and EN14683<sup>1</sup>
- ✓ Provide the **highest recognised level of protection from splash, offering ASTM Level 3 (160mmHg) protection**<sup>1</sup>
- ✓ NIOSH approved and **EXCEED NIOSH standards for filtration and breathability**<sup>2</sup>, designed to filter out at least 95% of particles at 0.3 microns
- ✓ **Entered in ARTG 351812** as both Medical Respirators and Surgical Masks
- ✓ **FDA cleared** Class 2 medical device
- ✓ **Chemotherapy protection.** Tested against 7 cytotoxic drugs<sup>3</sup>

# Mask Dispenser Promotion

**BUY 1 AND  
GET 1 FREE**

For a limited time only,  
while stocks last

**CLAIM**

**2 FOR 1  
SALE**




**TO CLAIM THIS OFFER, EMAIL PROOF OF PURCHASE TO  
[Communications\\_ANZ@hyh.com](mailto:Communications_ANZ@hyh.com).**

Terms and Conditions: 1. Offer available to customers in Australia and New Zealand only. 2. Buy 1 get 1 free applies to the same size dispenser only. 3. Applies to product codes 4966, 4967 and 4968 only. 4. To claim this offer, email proof of purchase to [Communications\\_ANZ@hyh.com](mailto:Communications_ANZ@hyh.com) during the promotional period. 5. For a limited time only. 6. Limits per customer apply. 7. While stocks last. 8. O&M Halyard Australia shall not be liable for any injury, loss or damage of any nature whatsoever (including but not limited to indirect or consequential losses) which is suffered or sustained as a result of or in connection with receiving, taking or using the products except for any liability which cannot be excluded by law. 9. If for any reason this offer is not capable of running as planned O&M Halyard Australia reserves the right in its sole discretion to modify the terms of offer including the products offered and these conditions of offer. 10. All other standard T&C's apply.

## More information

Contact your O&M Halyard Account Manager for more information about selecting an N95 Respirator that offers you the right level of protection.

1. IPS Job # 46161 and 72843
2. NIOSH Test Report TN-22659
3. CLAIM-07666. Code 46727 tested against (Carmustine, Cyclophosphamide, Doxorubican Hydrochloride, Etoposide, Fluorouracil, Paclitaxel, ThioTepa, with Carmustine breakthrough of 0.1 µg/cm<sup>2</sup>/min detected at 25.7 minutes);



O&M Halyard Australia

Level 1, Suite 1, Triniti 3 | 39 Delhi Rd | North Ryde NSW 2113

AU: 1800 664 227

NZ: 0800 18 17 16

[www.halyardhealth.com.au](http://www.halyardhealth.com.au)

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2021 O&M Halyard International. All Rights Reserved.

View our [Privacy Policy](#). [Unsubscribe](#).