




CONFIDENCE  
NOW COMES IN  
**POUCHES**



## **NEW HALYARD\*** Sterilization Pouches

Confidence in performance from  
the trusted brand in sterilisation wrap



Introducing the **NEW** range of **HALYARD\*** Sterilization Pouches

CONFIDENCE IN PERFORMANCE

Medical packaging serves as a sterile barrier system for instruments/equipment from sterilisation through to final use.

With the vast array of equipment you manage, in all their shapes and sizes, the last thing you need is a poor performing sterile barrier system that punctures and fails prior to presentation in theatres.

From the brand you know and trust in sterilisation wrap, Halyard now offers that same confidence in performance through our extensive range of reliable, high quality HALYARD\* Sterilization Pouches to meet your sterility needs.

With superior tear strength, puncture resistance, seal strength and an exceptional clean peel, the new HALYARD\* Sterilization Pouches provides confidence in performance you can see and feel.



Extensive range of sizes; heat and self-seal; flat pack and rolls; and gusseted packages.

From the brand you know and trust in sterilisation wrap



**Quick check of seal integrity** and performance with colour changing film<sup>1,2</sup>



**Clear-cut contrasting** process indicators where colour changes for:

STEAM	GREEN	PURPLE
EO	YELLOW	BROWN
H <sub>2</sub> O <sub>2</sub>	PINK	BLUE



**Puncture and abrasion resistance** with the specially formulated HALYARD\* multilayer transparent film<sup>1</sup>



**Excellent strength** under pressure with triple rail heat seals<sup>1</sup>



**Exceptional clean and fibre-free peel** from a specially designed seal between the multi-layer laminate and premium grade paper or Tyvek® backing<sup>1,3</sup>



**Compliance** to international norms and standards as well as AS/NZS 4187<sup>4</sup>





**Easy aseptic presentation** and opening with a top thumb notch<sup>1</sup>



Reduce collection of dust and **reduce the risk of contamination** upon opening with heat-sealed corners<sup>1</sup>



**Reduced risk of ink contamination** with all chemical indicators and printing external to the packing area<sup>1</sup>



**Outstanding sterilant penetration** and maintenance of bacterial barrier properties with premium 68gsm medical grade paper for steam & EO or Tyvek® 1073B material for H<sub>2</sub>O<sub>2</sub> gas plasma<sup>1,3,5</sup>

# Switch to HALYARD\* Sterilization Pouches for confidence in performance from a brand you know and trust.

For product samples or more  
information, please contact  
your local O&M Halyard  
representative or contact  
Customer Service on  
**1800 664 227** or email  
**CustomerService\_ANZ@hyh.com**

Product Code	Size	Compatible Sterilisation Modality	Packaging
--------------	------	-----------------------------------	-----------

## HEAT SEAL POUCHES

### HALYARD\* Heat-Seal Paper/Film Sterilization Pouches

48516	50x250mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48517	75x150mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48518	75x200mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48519	75x250mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48496	100x200mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48497	100x250mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48498	100x270mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48499	100x300mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48506	150x300mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48507	150x380mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48508	150x400mm	STEAM / EO	250 units/box; 8 boxes/case; 2000 units/case
48510	205x400mm	STEAM / EO	250 units/box; 4 boxes/case; 1000 units/case
48512	250x380mm	STEAM / EO	250 units/box; 4 boxes/case; 1000 units/case
48513	250x500mm	STEAM / EO	250 units/box; 2 boxes/case; 500 units/case

### HALYARD\* Heat-Seal Tyvek®/Film Sterilization Pouches

48532	75x210mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case
48527	100x260mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case
48528	150x320mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case
48529	200x400mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case
48530	200x500mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case
48531	200x650mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	250 units/box; 8 boxes/case; 2000 units/case

### HALYARD\* Heat Seal Paper/Film Sterilization Reels

48525	50mmx200m	STEAM / EO	2 reels/case
48526	75mmx200m	STEAM / EO	2 reels/case
48521	100mmx200m	STEAM / EO	1 reel/case
48522	150mmx200m	STEAM / EO	1 reel/case
48523	200mmx200m	STEAM / EO	1 reel/case
48524	250mmx200m	STEAM / EO	1 reel/case

### HALYARD\* Heat-Seal Tyvek®/Film Sterilization Reels

48534	150mm x 100m	PLASMA (H <sub>2</sub> O <sub>2</sub> )	1 reel/box; 8 boxes/case; 8 reels/case
48535	200mm x 100m	PLASMA (H <sub>2</sub> O <sub>2</sub> )	1 reel/box; 4 boxes/case; 4 reels/case
48536	250mm x 100m	PLASMA (H <sub>2</sub> O <sub>2</sub> )	1 reel/box; 4 boxes/case; 4 reels/case

### HALYARD\* Heat-Seal Paper/Film Gusseted Sterilization Pouches

48493	200 x 55 x 400mm	STEAM / EO	100 units/box; 10 boxes/case; 1000 units/case
48494	150mm x 50mm x 100m	STEAM / EO	2 reels/case

## SELF SEAL POUCHES

### HALYARD\* Self-Seal Paper/Film Sterilization Pouches

48548	90 x 130mm	STEAM / EO	200 units/box; 10 boxes/case; 2000 units/case
48549	90 x 200mm	STEAM / EO	200 units/box; 10 boxes/case; 2000 units/case
48550	90 x 250mm	STEAM / EO	200 units/box; 10 boxes/case; 2000 units/case
48541	200 x 350mm	STEAM / EO	200 units/box; 5 boxes/case; 1000 units/case
48544	250 x 400mm	STEAM / EO	200 units/box; 5 boxes/case; 1000 units/case
48545	300 x 450mm	STEAM / EO	200 units/box; 5 boxes/case; 1000 units/case

### HALYARD\* Self-Seal Tyvek®/Film Sterilization Pouches

48552	150 x 320mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	200 units/box; 10 boxes/case; 2000 units/case
48554	75 x 200mm	PLASMA (H <sub>2</sub> O <sub>2</sub> )	200 units/box; 10 boxes/case; 2000 units/case

Additional sizes available on request



Infection Prevention Product Solutions  
Halyard\* Education Foundation  
Tools & Best Practice  
Commitment to Excellence and Sustainability  
Experienced Customer Support  
Knowledgeable Sales Force  
Innovation & Technology

#### References

1. Data on file. Claim-07729
2. Paper/Film pouches only
3. Billerudkorsnäs Technical Data Sheet
4. ISO 11607, EN868, ANSI/AAMI ST79, AS/NZS 4187
5. Tyvek 1073B Medical Device Packaging <https://www.dupont.com/products/tyvek-1073b-medical-device-packaging.html>

