



SURGICAL PRODUCT
CATALOGUE

AUSTRALIA &
NEW ZEALAND

QUICK LINKS

SURGICAL GOWNS

AERO CHROME*	page 50
AERO BLUE*	page 51
Surgical Gowns	page 52

SURGICAL DRAPES & PACKS

Cardiac Catheterisation / Radiology Drapes	page 53
Cardiovascular Drapes	page 55
ENT Drapes & Packs	page 56
Equipment Covers	page 58
General Drapes & Packs:	
Accessories	page 59
Basic Drapes	page 60
Breast / Chest Drapes	page 62
Fenestrated Drapes	page 63
Thyroid Drapes	page 63
Laparoscopy / Laparotomy Drapes	page 64
Minor Procedure Drapes	page 65
Universal Packs	page 66
Universal Drapes	page 69
Neurosurgical Drape	page 71
Obstetrics and Gynaecology Drapes & Packs:	
C-Section Drapes	page 72
Laparoscopy	page 73
Lithotomy Packs	page 74
Lithotomy Drapes	page 75
Under Buttocks Drapes	page 76

Ophthalmic Drapes & Packs	page 77
--------------------------------------	---------

Orthopaedic Drapes & Packs:

Stockinettes	page 78
Ortho Pack	page 78
Extremity Drapes	page 79
Hand Drapes	page 81
Hip Drapes	page 82
Shoulder Drapes	page 83
U Drapes	page 85
Bar Drapes	page 87
Knee Arthroscopy Drapes	page 88

Urology Drapes & Packs	page 89
-----------------------------------	---------

STANDARD PACKS

ENT Pack	page 90
Universal Pack	page 90
<u>Cysto-Gynae Pack</u>	page 91
Caesarean Pack	page 91
Lower Extremity Pack	page 92

CUSTOM PROCEDURE PACKS (CPP)

	page 93
--	---------

SURGICAL GOWNS



AERO CHROME* BREATHABLE PERFORMANCE SURGICAL GOWN

- Clearly Identifiable Silver Colour
- Soft, Light and Breathable Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 4 / EN 13795 High Performance

Code	Description	Packaging
48578	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Small	Items/Case = 30
48579	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large	Items/Case = 30
48392	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
48583	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
48399	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
48580	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
48393	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
48584	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
48400	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
48581	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
48394	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
48585	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
48401	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
48582	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 28

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products..

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

SURGICAL GOWNS



AERO BLUE* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24


All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

SURGICAL GOWNS



SURGICAL GOWNS

- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm


Code	Description	Packaging
90003	Surgical Gown, Sterile, Small–Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large	Items/Case = 640
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large	Items/Case = 640
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
Contents:		
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIAC CATHETERISATION / RADIOLOGY DRAPES



PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm


Code	Description	Packaging
89368	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12



BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm

- CONTROL* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration

Code	Description	Packaging
89471	Brachial Angiography Drape, Sterile	Items/Case = 50



FEMORAL ANGIOGRAPHY / IJ ACCESS COMBINATION DRAPES, 200cm x 360cm

- Two 17cm Trapezoidal Fenestrations
- Two 9cm Adhesive Circular Fenestrations
- Two 58cm x 360cm Clear Windows

Code	Description	Packaging
29225	Femoral Angiography / IJ Access Combination Drape with Windows, Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIAC CATHETERISATION / RADIOLOGY DRAPE



RADIAL FEMORAL ANGIOGRAPHY DRAPE, 244cm x 381cm

- Two 7cm Adhesive Circular Fenestrations
- Two 4cm x 6cm Oval Fenestrations
- Two 61cm x 381cm Clear Windows
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44976	Radial/Femoral Angiography Drape, Sterile	Items/Case = 6



EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
89749	Extra-Long Femoral Angiography Drape, Sterile	Items/Case = 16
79749	Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIOVASCULAR DRAPES



CVARTS* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89452	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7



CVARTS* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89453	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18




CVARTS* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89454	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5
Contents:		
1-CVARTS* Top Drape, 254cm x 152cm		
1-CVARTS* Bottom Drape, 193cm x 152cm		
2-CVARTS* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop		


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ENT DRAPES & PACKS




ENT PACK III

Code	Description	Packaging	
88223	ENT Pack III	Items/Case = 6	
Sterile Component Code	Non-Sterile Component Code	Contents	
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm	
89431	79430	1-Mayo Stand Cover 80cm x 145cm	
89352	79351	1-Suture Bag	
		1-Head Drape 91cm x 127cm:	
		• Turban	
		• Positioning Tape	
		1-ENT Split Drape 193cm x 315cm:	
		• with 7 x 71cm Split	
		• Positioning Tape	
		• CONTROL* Fabric Reinforcement	
		• Tube Holders	



EENT PACK IV

Code	Description	Packaging	
88719	EENT Pack IV	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
89701	79720	1-Surgical Gown with Set-in Sleeves, In Overwrap, with Towel, Large	
89601	79601	1-Surgical Gown, with Set-in Sleeves, Large	
89438	79438	1-Absorbent Towel, 38cm x 56cm	
		1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
		1-Suture Bag	
		1-Head Drape, 104cm x 81cm	
		1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split	



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad

Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ENT DRAPE & PACKS



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm


Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38




DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm from the top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

EQUIPMENT COVERS



MAYO STAND COVERS

Code	Description	Packaging
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44



BACK TABLE COVERS

Code	Description	Packaging
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20

OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
89621	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS

ACCESSORIES




ABSORBENT TOWELS, 38cm x 56cm

Code	Description	Packaging
89701	Absorbent Towels, White, Sterile, Pack of 2	Items/Case = 100
79720	Absorbent Towels, White, Non-Sterile	Items/Case = 1000



TAPE STRIPS, 10cm x 55cm

Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300




UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS

BASIC DRAPES




BASICS* UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89678	Utility Drape, Sterile	Items/Case = 30
79782	Utility Drape, Non-Sterile	Items/Case = 90




HALF DRAPE, 147cm x 101cm

Code	Description	Packaging
47614	Half Drape, Sterile	Items/Case = 54
47615	Half Drape, Non-Sterile	Items/Case = 1000



MEDIUM DRAPE, 185cm x 101cm

Code	Description	Packaging
47617	Medium Drape, Sterile	Items/Case = 60
47619	Medium Drape, Non-Sterile	Items/Case = 800



MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS



LARGE DRAPE, 139cm x 193cm


Code	Description	Packaging
47622	Large Drape, Sterile	Items/Case = 36
47624	Large Drape, Non-Sterile	Items/Case = 600



THREE-QUARTER DRAPE, 139cm x 193cm


- CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
47628	Three-Quarter Drape, Sterile	Items/Case = 20



X-LARGE DRAPE, 193cm x 243cm


Code	Description	Packaging
47630	X-Large Drape, Sterile	Items/Case = 24
47631	X-Large Drape, Non-Sterile	Items/Case = 24



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


BREAST / CHEST DRAPES



CHEST U-DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



BREAST DRAPE, 254cm x 193cm x 315cm


- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


FENESTRATED DRAPES



EPIDURAL DRAPE WITH TAPE, 50cm x 75cm

- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
89636	HALYARD BASICS* Epidural Drape with Tape, Sterile	Items/Case = 40



VARIABLE FENESTRATED DRAPE, 75cm x 100cm

- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
89680	HALYARD BASICS* Variable Fenestrated Drape, Sterile	Items/Case = 30



CV LINE DRAPE, 160cm x 240cm

- CONTROL* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code	Description	Packaging
89835	CV Line Drape, Sterile	Items/Case = 10
79835	CV Line Drape, Non-Sterile	Items/Case = 14



THYROID DRAPE, 193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders

Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers


Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


MINOR PROCEDURE DRAPES



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Rectangular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Circular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Rectangular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


UNIVERSAL PACKS

UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



Code	Description	Packaging
88759	Universal Pack – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 38cm x 56cm 1-Suture Bag
34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm

UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES



Code	Description	Packaging
88769	Universal Pack – Minor	Items/Case = 14
Sterile Component Code	Non-Sterile Component Code	Contents
89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
	79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
	79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
	79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS




UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
88774	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 2-Absorbent Towels, 38cm x 56cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS



U-BAR* PACK

Code	Description	Packaging	
88334	U-BAR* Pack	Items/Case = 16	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm	



U-BAR* PACK – IV

Code	Description	Packaging	
88773	U-BAR* Pack – IV	Items/Case = 10	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


UNIVERSAL DRAPES



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS



TOP DRAPE, 276cm x 150cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70




SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

NEUROSURGICAL DRAPE



CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm

- CONTROL* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port


Code	Description	Packaging
89062	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


C-SECTION DRAPES



C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89528	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89529	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LAPAROSCOPY

LAPAROSCOPY / PELVISCOPY PACK I




Code Description

88222 Laparoscopy/Pelviscopy Pack

Items/Case = 6

Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm
		1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm
		• 18cm x 23cm Trapezoidal Abdominal Fenestration
		• 8cm x 15cm Perineal Fenestration with Flap



LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers

Code Description

89219 Laparoscopy/Pelviscopy Drape Drape with Pouch, Sterile

Packaging

Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY


DRAPES & PACKS

LITHOTOMY PACKS

LITHOTOMY PACK II



Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
		1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
47614	79101	1-Half Drape, 147cm x 101cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LITHOTOMY DRAPES



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY


DRAPES & PACKS

UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30




UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7cm Triangular Pouch

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


OPHTHALMIC DRAPE & PACKS



BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



OCUARTS* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders

Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS



STOCKINETTES

Code	Description	Packaging
89811	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
48759	Stockinette, Non-Sterile, Medium, 23cm x 112cm	Items/Case = 90
89821	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
48760	Stockinette, Non-Sterile, Large, 30cm x 122cm	Items/Case = 56



ORTHO PACK I – MAJOR

Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
44682	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
89331	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split
		1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
89381	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


EXTREMITY DRAPES



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84




ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPE & PACKS



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE, 193cm x 304cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders

Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBOARD COVERS, 295cm x 223cm x 366cm


- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


HAND DRAPES



HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
77450	Hand Drape, Sterile	Items/Case = 11
13227	Hand Drape, Non-Sterile	Items/Case = 108



ORTHOARTS* HAND DRAPE, 193cm x 381cm


- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


HIP DRAPES



ORTHOARTS* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89345	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



HIP DRAPE WITH CLEAR LEG POCKETS, 294cm 223cm x 345cm


- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
89349	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders


Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric

Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


U DRAPES



U-DRAPE, 193cm x 305cm

- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190



LARGE U-DRAPE, 223cm x 317cm

- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15




ORTHOARTS* U-DRAPE, 223cm x 274cm

- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* BODY SPLIT DRAPE, 223cm x 294cm

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12



SPLIT DRAPE, 152cm x 178cm

- Impervious
- Polyethylene Material
- 10cm x 53cm Split

Code	Description	Packaging
89331	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



LARGE SPLIT DRAPE, 152cm x 205cm

- Impervious
- Polyethylene Material
- 11cm x 78cm Split

Code	Description	Packaging
89333	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


BAR DRAPES



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28



VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14



ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
89093	Isolation Screen, Sterile	Items/Case = 24

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


KNEE ARTHROSCOPY DRAPES



ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7




ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

UROLOGY DRAPES & PACKS



T.U.R. PACK

Code	Description	Packaging
77885	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component Code	Contents
		1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Absorbent Towel, 38cm x 50cm
76885	1-T.U.R. Drape, 289cm x 175cm x 170cm:	
		<ul style="list-style-type: none"> • Attached Leggings • Fluid Collection Pouch with Screen, Suction Port and Tie • 54cm x 75cm Triangular Pouch • Finger Cot • 7cm Fenestration • 25cm x 20cm Fenestration



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
	Sterile Component Code	Non-Sterile Component Code
		Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm with 7cm x 15cm Perineal Fenestration and 23cm x 23cm Screen

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

STANDARD PACKS

Please contact your O&M Halyard Account Manager for pack information.

ENT PACK

Description	Packaging
ENT Pack	Items/Case = 1
<hr/>	
Contents	
<hr/>	
1-Back Table Cover, Unfolded, 152cm x 229cm	
1-Mayo Stand Cover, 58cm x 137cm	
1-Drape, Head Turban, 127cm x 91cm	
1-Drape, EENT, Split, 193cm x 315cm	
2-Blade, Surgical, Carbon Steel, 15	
5-Gauze, 8 Ply, 10cm x 10cm, XRD	
1-Cover, Light Handle, Flexible, Green	
2-Syringe, 10ml, L/L	
2-Syringe, 20ml, L/L	
2-Needle, 18G x 38mm	
2-Needle, 23G x 25mm	
1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
1-Tubing, Suction, Soft, 4m	
1-Handle, Yankaeur, Bulb Tip, Non-Vented	
3-Gallipot, 150ml	
2-Dish, Kidney, 700ml, Clear	
1-Bowl, 250ml, Clear	
1-Bowl, 1000ml	
1-Tray, Sharps Transfer, Yellow	

UNIVERSAL PACK

Description	Packaging
Universal Pack	Items/Case = 1
<hr/>	
Contents	
<hr/>	
1-Back Table Cover, Unfolded, 152cm x 229cm	
1-Mayo Stand Cover, 58cm x 137cm	
1-Drape, Top, 273 x 127cm	
1-Drape, Bottom, 193cm x 195cm	
2-Drape, Side, 195cm x 105cm	
1-Blade, Surgical, Carbon Steel, 15	
1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
5-Gauze, 8 Ply, 10cm x 10cm, XRD	
5-Sponge, Abdominal, XRD, 30cm x 30cm	
1-Cover, Light Handle, Flexible, Green	
1-Marker, Skin, Standard Tip	
1-Ruler, White, 15cm	
1-Syringe, 20ml, L/L	
2-Needle, 23G x 25mm	
1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
5-Label, 125mm x 5cm	
2-Gallipot, 150ml	
2-Dish, Kidney, 700ml, Clear	
1-Bowl, 1000ml	
1-Tray, Sharps Transfer, Yellow	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

STANDARD PACKS

Please contact your O&M Halyard Account Manager for pack information.

CYSTO-GYNAE PACK

Description	Packaging
Cysto-Gynae Pack	Items/Case = 1
Contents	
1-Back Table Cover, Unfolded, 152cm x 229cm	
1-Drape, Lithotomy T	
2-Drape, Leggings, 78cm x 122cm	
5-Gauze, 8 Ply, 10cm x 10cm, XRD	
1-Pad, Peri Winged 10cm x 27cm	
2-Gallipot, 150ml	
1-Dish, Kidney, 700ml, Clear	
1-Dish, Kidney, 700ml, Yellow	
1-Bowl, 1000ml	

CAESAREAN PACK

Description	Packaging
Caesarean Pack	Items/Case = 1
Contents	
1-Back Table Cover, Unfolded, 152cm x 229cm	
1-Mayo Stand Cover, 58cm x 137cm	
1-Drape, C-Section, Fenestrated Incise	
1-Bag, Suture, with tape, 20cm x 25cm	
1-Blade, Surgical, Carbon Steel, 23	
1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
5-Gauze, 8 Ply, 10cm x 10cm, XRD	
5-Sponge, Abdominal, XRD, 30cm x 30cm	
1-Cover, Light Handle, Flexible, Green	
1-Syringe, 10ml, L/S	
1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
1-Tubing, Suction, 3.6m	
1-Handle, Yankauer, Bulb Tip, Non-Vented	
1-Catheter, Y-Suction, 10FG x 45cm	
2-Gallipot, 150ml	
2-Dish, Kidney, 700ml, Clear	
1-Bowl, 1000ml	
1-Tray, Sharps Transfer, Yellow	
1-Blanket, Baby, 75cm x 100cm	
1-Pad, Peri Winged 10cm x 27cm	

NOTE: Please see package inserts for directions for use for each of these products..
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

STANDARD PACKS

Please contact your O&M Halyard Account Manager for pack information.

LOWER EXTREMITY PACK

Description	Packaging
Lower Extremity Pack	Items/Case = 1
Contents	
1-Back Table Cover, Unfolded, 152cm x 229cm	
1-Drape, Extremity, Extra Long, 289cm x 321cm	
1-Drape, Split, 152cm x 178cm	
1-Stockinette, Impervious, Large, 30cm x 122cm	
1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
5-Gauze, 8 Ply, 10cm x 10cm, XRD	
5-Sponge, Abdominal, XRD, 30cm x 30cm	
1-Bandage, 15cm x 2.7m, Medium Weight	
1-Cover, Light Handle, Flexible, Green	
1-Marker, Skin, Standard Tip	
1-Ruler, White, 15cm	
5-Label, 125mm x 500mm, Blank	
1-Syringe, 20ml, L/L	
1-Needle, 23G x 38mm	
1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
2-Gallipot, 150ml	
2-Dish, Kidney, 700ml, Clear	
1-Bowl, 1000ml	
1-Tray, Sharps Transfer, Yellow	

NOTE: Please see package inserts for directions for use for each of these products..
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

CUSTOM PROCEDURE PACKS (CPP)

THE FUTURE OF CUSTOM PROCEDURE PACKS FROM THE NAME YOU TRUST

For over 30 years HALYARD has been providing high quality surgical drapes and gowns, designed to deliver protection, comfort and optimal patient outcomes.



Leveraging Halyard's global expertise in supply chain, component sourcing, utilisation reviews and cost efficiencies, we are now extending our offering to provide ANZ clinicians with HALYARD branded Custom Procedure Packs (CPPs).

Based on years of working together with clinicians globally we know that some of the challenges you face in Operating Theatres include:

- Operational inefficiencies
- Lack of visibility to costs and inventories
- Suboptimal utilisation of Operating Theatre supplies
- Procedural supply waste
- Need for improved throughput
- Staff satisfaction and retention

That's where HALYARD comes in. When you make the decision to switch to HALYARD CPPs, our clinical experts work with you to ensure EVERY component is used EVERY time. It's part of our PACT pledge and the HALYARD CPP Solution.

NOTE: Please see package inserts for directions for use for each of these products..
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.