



PATIENT CARE  
ACCESSORIES  
CATALOGUE

AUSTRALIA &  
NEW ZEALAND



# QUICK LINKS

## COLD THERAPY

### Ice Packs

HALYARD* Bilateral Facial Ice Pack	page 3
HALYARD* EyeCare (EENT) Ice Pack	page 3
HALYARD* STAY-DRY Ice Pack	page 3

## I.V. SUPPORT

### I.V. Arterial And Wrist Supports

HALYARD* Neonatal & Infant/Child I.V. Support	page 4
HALYARD* HAND-AID* Arterial Wrist Support	page 4
HALYARD* HAND-AID* Wrist Support	page 4
HALYARD* Antecubital I.V. Support	page 4

## O.R. ACCESSORIES

### Support Straps

HALYARD* Stirrup Strap	page 5
HALYARD* TS-30* Personal Utility Strap	page 5
HALYARD* Tourniquet Cover	page 5

## PATIENT LINEN

### Single Use Linen

HALYARD* Stretcher Sheet	page 6
HALYARD* Draw Sheet	page 6
HALYARD* Pillowcase	page 6

# COLD THERAPY

## ICE PACKS



### HALYARD\* BILATERAL FACIAL ICE PACK

- 2 ice packs bonded together to evenly deliver cold therapy to both sides of the jaw
- Can be used for recovery after maxillofacial and dental procedures
- Single-patient use

Code	Description	Size	Packaging
33101	Bilateral Facial Ice Pack	13cm x 30cm	12 eaches/box 2 boxes/case 24 eaches/case



### HALYARD\* EYECARE (EENT) ICE PACK

- Designed for easy application after plastic surgery, ophthalmic and sinus procedures
- Adjustable strap ensures the pack sits comfortably and securely on the face
- Single-patient use

Code	Description	Size	Packaging
33150	EyeCare (EENT) Ice Pack, 1 Strap	11.4cm x 25cm	15 eaches/box 2 boxes/case 30 eaches/case



### HALYARD\* STAY-DRY ICE PACK

- For general use
- Delivers up to 2.5 hours of consistent, dry cold therapy solutions
- STAY-DRY construction features 3 layers which help to reduce wetness due to condensation
- 4 ties for versatile application on shoulders, knees, elbows, etc
- Single-patient use

Code	Description	Size	Packaging
33495	STAY-DRY Ice Pack, Small, 4 Ties	13cm x 30cm	25 eaches/box 2 boxes/case 50 eaches/case
33500	STAY-DRY Ice Pack, Large, 4 Ties	16.5cm x 36cm	25 eaches/box 2 boxes/case 50 eaches/case

# I.V. SUPPORT

## I.V. ARTERIAL AND WRIST SUPPORTS



### HALYARD\* NEONATAL & INFANT/CHILD I.V. SUPPORT

- Padded for comfort
- Aluminium core easily bends to conform to various angles of small wrists and feet
- Single-patient use

Code	Description	Size	Packaging
29965	Neonatal I.V. Support	9cm x 4cm, 13cm straps	30 eaches/case
29967	Infant/Child I.V. Support	16cm x 5cm, 23cm straps	30 eaches/case



### HALYARD\* HAND-AID\* ARTERIAL WRIST SUPPORT

- Uniquely anatomically designed to position the wrist for accurate and efficient application of arterial lines
- The instant support of the foam block saves setup time
- Single-patient use

Code	Description	Size	Packaging
29990	Hand-Aid* Arterial Wrist Support, Paediatric	16cm x 7cm, 23cm straps	20 eaches/case
29980	Hand-Aid* Arterial Wrist Support, Adult	23cm x 9cm, 37cm straps	20 eaches/case
29985	Hand-Aid* Arterial Wrist Support, Replacement Straps & Starter Block, Compatible with 29980		50 eaches/case



### HALYARD\* HAND-AID\* WRIST SUPPORT

- Anatomically designed to position the wrist for accurate and efficient cannulisation in dorsal hand veins
- Single-patient use

Code	Description	Size	Packaging
29970	Hand-Aid* I.V. Wrist Support, Paediatric	14cm x 5cm, 23cm straps	20 eaches/case
29900	Hand-Aid* I.V. Wrist Support, Adult	16cm x 7cm, 23cm straps	20 eaches/case



### HALYARD\* ANTECUBITAL I.V. SUPPORT

- The unique and comfortable way to immobilise and support during antecubital I.V. therapy
- Single-patient use

Code	Description	Size	Packaging
29920	Antecubital I.V. Support	32cm x 11cm, 48 cm straps	10 eaches/case

# O.R. ACCESSORIES

## SUPPORT STRAPS



### HALYARD\* STIRRUP STRAP

- Designed to ensure patient comfort in the lithotomy position
- Made of soft fabric and featuring a layer of padding
- Comes with slip ring for easy adjustment
- Single-patient use

Code	Description	Size	Packaging
52712	Personal Stirrup Strap with Slip Ring	48cm x 9cm	10 eaches/pack 4 packs/case 40 eaches/case



### HALYARD\* TS-30\* PERSONAL UTILITY STRAP

- Designed to secure and position the patient during a surgical procedure
- Hook and loop closure for easy adjustment of the strap
- Soft fabric for patient comfort
- Single-patient use

Code	Description	Size	Packaging
55299	TS-30* Personal Utility Strap	76cm x 8cm	12 eaches/pack 2 packs/case 24 eaches/case



### HALYARD\* TOURNIQUET COVER

- Designed to cover a reusable, pneumatic tourniquet cuff to protect it from skin prep and other procedural fluids
- Soft padding helps prevent pinching and wrinkling of patients' skin
- Comes in two sizes – small and large
- Single-patient use

Code	Description	Size	Packaging
53600	O.R. Tourniquet Cover, Small	53.5cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case
53605	O.R. Tourniquet Cover, Large	78cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case

# PATIENT LINEN



## HALYARD\* SINGLE USE LINEN

<b>Code</b>	<b>Description</b>	<b>Size</b>	<b>Packaging</b>
67783	Stretcher Sheet, Blue	101.5cm x 213cm (40" x 84"), Bulk	80 eaches/case
67843	Draw Sheet, Poly-Coated*, Blue <small>*Coated with Polyethylene Film</small>	101.5cm x 183cm (40" x 72"), Bulk	80 eaches/case
67803	Pillowcase, Blue	50.5cm x 73.5cm (20in x 29in), Bulk	200 eaches/case



Infection Prevention Product Solutions  
Halyard\* Education Foundation  
Tools & Best Practice  
Commitment to Excellence and Sustainability  
Experienced Customer Support  
Knowledgeable Sales Force  
Innovation & Technology



**O&M Halyard Australia Pty. Limited**  
Customer Care:  
Phone: 1800 664 227  
Email: CustomerService\_ANZ@hyh.com

**[www.halyardhealth.com.au](http://www.halyardhealth.com.au)**

**[www.halyardhealth.com.au](http://www.halyardhealth.com.au)**



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2022 O&M Halyard, Inc.  
All Rights Reserved. 22:06:WFHY947v2