



SURGICAL PRODUCT CATALOGUE

AUSTRALIA &
NEW ZEALAND

QUICK LINKS

SURGICAL GOWNS

AERO CHROME*	page 3
AERO BLUE*	page 4
Surgical Gowns	page 5

SURGICAL DRAPES & PACKS

Cardiac Catheterisation / Radiology Drapes	page 6	Ophthalmic Drapes & Packs	page 35
Cardiovascular Drapes	page 8	Orthopaedic Drapes & Packs:	
ENT Drapes & Packs	page 9	Ortho Packs	page 36
Equipment Covers	page 12	Extremity Drapes	page 38
General Drapes & Packs:		Hand Drapes	page 40
Accessories	page 13	Hip Drapes	page 41
Basic Packs	page 14	Knee Arthroscopy Drapes	page 42
Basic Drapes	page 15	Shoulder Drapes	page 44
Breast / Chest Drapes	page 17	U & Bar Drapes	page 46
Fenestrated Drapes	page 18	Urology Drapes & Packs	page 50
Laparoscopy / Laparotomy Drapes	page 19		
Minor Procedure Drapes	page 21		
Thyroid Drapes	page 22		
Universal Packs	page 23		
Universal Drapes	page 26		
Neurosurgical Drape	page 28		
Obstetrics and Gynaecology Drapes & Packs:			
C-Section Drapes	page 29		
Laparoscopy Packs	page 30		
Laparoscopy Drapes	page 31		
Lithotomy Packs	page 32		
Lithotomy Drapes	page 33		
Under Buttocks Drapes	page 34		
		STANDARD PACKS	
		ENT Pack	page 51
		Universal Pack	page 51
		Cysto-Gynae Pack	page 52
		Caesarean Pack	page 52
		Lower Extremity Pack	page 53
		CUSTOM PROCEDURE PACKS (CPP)	page 54
		INDEX	page 55

SURGICAL GOWNS



AERO CHROME* BREATHABLE PERFORMANCE SURGICAL GOWN

- Soft Knitted V-neck Collar
- Clearly Identifiable Silver Colour
- Soft, Light and Breathable Fabric
- Raglan Sleeves
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 4 / EN 13795 High Performance

Code	Description	Packaging
44680	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Small	Items/Case = 30
44661	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
44681	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large	Items/Case = 30
44662	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
44682	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
44664	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
44683	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
44666	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
44684	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24
44685	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
44667	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
44686	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
44668	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
44687	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
44669	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.. Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

SURGICAL GOWNS



AERO BLUE* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
43941	AERO BLUE* Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
43623	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24


All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

SURGICAL GOWNS



SURGICAL GOWNS

- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- EN 13795 Standard Performance

Code	Description	Packaging
90003	Surgical Gown, Sterile, Small–Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large	Items/Case = 640
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large	Items/Case = 640
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
70392	Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 60
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
	Contents:	
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	


All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIAC CATHETERISATION / RADIOLOGY DRAPES



PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm


Code	Description	Packaging
89368	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12



BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm

- CONTROL* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration

Code	Description	Packaging
89471	Brachial Angiography Drape, Sterile	Items/Case = 50



RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240cm x 380cm

- Two 63cm x 350cm Windows
- Two 12cm x 7cm Oval Radial Fenestrations
- Two 12cm Round Femoral Fenestrations

Code	Description	Packaging
89583	HALYARD BASICS* Radial Femoral Angiography Drape, Sterile	Items/Case = 12

RADIAL FEMORAL ANGIOGRAPHY DRAPE, 244cm x 381cm

- Two 7cm Adhesive Circular Fenestrations
- Two 4cm x 6cm Oval Fenestrations
- Two 61cm x 381cm Clear Windows
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44976	Radial/Femoral Angiography Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIAC CATHETERISATION / RADIOLOGY DRAPES



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOW, 208cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- One 68cm x 315cm Windows
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89709	Extra Long Femoral Angiography Drape with Window, Sterile	Items/Case = 14
79709	Extra Long Femoral Angiography Drape with Window, Non-Sterile	Items/Case = 16



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOWS, 221cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- Two 68cm x 343cm Windows
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89711	Extra Long Femoral Angiography Drape with Dual Windows, Sterile	Items/Case = 16
13296	Extra Long Femoral Angiography Drape with Dual Windows, Non-Sterile	Items/Case = 180



EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
89749	Extra-Long Femoral Angiography Drape, Sterile	Items/Case = 16
79749	Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


CARDIOVASCULAR DRAPES



PERINEAL DRAPE, 187cm x 241cm

- Split 44cm x 95cm
- CONTROL* Plus Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
76356	Perineal Drape, Non-Sterile	Items/Case = 17



CVARTS* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89454	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5
Contents:		
1-CVARTS* Top Drape, 254cm x 152cm		
1-CVARTS* Bottom Drape, 193cm x 152cm		
2-CVARTS* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop		



CVARTS* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89452	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7



CVARTS* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89453	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric





Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ENT DRAPES & PACKS




EENT PACK I

Code	Description	Packaging
88711	EENT Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown in Overwrap, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Absorbent Towel, 38cm x 56cm
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split

EENT PACK II

Code	Description	Packaging
88716	EENT Pack II	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown In Overwrap, Fabric Reinforced, with Towel, Large
95211	74911	1-ULTRA Surgical Gown, Fabric Reinforced, with Towel, Large
		1-Microsurgery Mayo Stand Cover, 58cm x 135cm
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



ENT PACK III

Code	Description	Packaging
88223	ENT Pack III	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Mayo Stand Cover 80cm x 145cm
		1-Suture Bag
89431	79430	1-Head Drape 91cm x 127cm:
		• Turban
		• Positioning Tape
89352	79351	1-ENT Split Drape 193cm x 315cm:
		• with 7 x 71cm Split
		• Positioning Tape
		• CONTROL* Fabric Reinforcement
		• Tube Holders

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

EENT DRAPES & PACKS



EENT PACK IV


Code	Description	Packaging
88719	EENT Pack IV	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Set-in Sleeves, In Overwrap, with Towel, Large
89701	79720	1-Surgical Gown, with Set-in Sleeves, Large
89601	79601	1-Absorbent Towel, 38cm x 56cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89438	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm

Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ENT DRAPES & PACKS



DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm from the top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15




EENT FENESTRATED DRAPE, 183cm x 229cm

- 6cm x 5cm Fenestration
- CONTROL Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89421	EENT Fenestrated Drape, Sterile	Items/Case = 18


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

EQUIPMENT COVERS



MAYO STAND COVERS

Code	Description	Packaging
88665	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Sterile	Items/Case = 80
78656	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Non-Sterile	Items/Case = 140
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



BACK TABLE COVERS

Code	Description	Packaging
88666	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Sterile	Items/Case = 45
78664	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Non-Sterile	Items/Case = 60
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300
37498	Back Table Cover, Two Tier Heavy Duty, 152cm, Sterile	Items/Case = 10

OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
89621	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS

ACCESSORIES




ABSORBENT TOWELS, 38cm x 56cm

Code	Description	Packaging
89701	Absorbent Towels, White, Sterile, Pack of 2	Items/Case = 100
79720	Absorbent Towels, White, Non-Sterile	Items/Case = 1000



TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300



DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
89945	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS

BASIC PACKS

BASIC PACK VI

Code	Description	Packaging
88161	Basic Pack VI	Items/Case = 7
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap with Towel, Large
95121	71321	1-ULTRA Surgical Gown, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
89731	79731	4-Utility Drapes, 66cm x 38cm, with Tape
		1-Suture Bag
89101	79101	1-Half Drape, 152cm x 112cm
89121	79121	1-Large Drape, 152cm x 193cm

BASIC PACK XV

Code	Description	Packaging
88178	Basic Pack XV	Items/Case = 15
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
89121	79121	1-Large Drape, 152cm x 193cm

SET-UP PACK I

Code	Description	Packaging
88171	Set-Up Pack I	Items/Case = 14
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-Ultra Surgical Gown In Overwrap with Towel, Large
95121	71321	1-ULTRA Surgical Gown, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS

BASIC DRAPES




UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89353	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




HALF DRAPE, 152cm x 101cm

Code	Description	Packaging
47614	Half Drape, Sterile	Items/Case = 54
47615	Half Drape, Non-Sterile	Items/Case = 1000




MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
89111	Medium Drape, Sterile	Items/Case = 50



MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50



MEDIUM DRAPE, 185cm x 101cm

Code	Description	Packaging
47619	Medium Drape, Non-Sterile	Items/Case = 800

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS



LARGE DRAPE, 139cm x 193cm


Code	Description	Packaging
47522	Large Drape, Sterile	Items/Case = 36
47624	Large Drape, Non-Sterile	Items/Case = 600



THREE-QUARTER DRAPE, 152cm x 193cm


- 39CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
47628	Three-Quarter Drape, Sterile	Items/Case = 20



X-LARGE DRAPE, 193cm x 243cm


Code	Description	Packaging
47630	X-Large Drape, Sterile	Items/Case = 20
47631	X-Large Drape, Non-Sterile	Items/Case = 24



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


BREAST / CHEST DRAPES



CHEST U-DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



BREAST DRAPE, 254cm x 193cm x 315cm


- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS

FENESTRATED DRAPES



EPIDURAL DRAPE WITH TAPE, 50cm x 75cm


- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape

Code Description

89636 HALYARD BASICS* Epidural Drape with Tape, Sterile

Packaging

Items/Case = 40



VARIABLE FENESTRATED DRAPE, 75cm x 100cm


- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape

Code Description

89680 HALYARD BASICS* Variable Fenestrated Drape, Sterile

Packaging

Items/Case = 30



CV LINE DRAPE, 160cm x 240cm

- CONTROL* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code Description


89835 CV Line Drape, Sterile

79835 CV Line Drape, Non-Sterile

Packaging

Items/Case = 10

Items/Case = 14



CV FULL BODY DRAPE, 172cm x 279cm

- CONTROL* Plus Fabric Reinforcement
- 10cm Circular Adhesive

Code Description

31445 CV Full Body Drape, Non-Sterile

Packaging

Items/Case = 27

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers


Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



LAPAROSCOPIC DRAPE, 264cm x 89cm x 315cm

- 31cm x 33cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Clear Film Pouch
- Tube Holders
- Armboard Covers
- Under Leg Wings

Code	Description	Packaging
29240	Laparoscopic Drape, Sterile	Items/Case = 15
29241	Laparoscopic Drape, Non-Sterile	Items/Case = 19



HALYARD BASICS* LAPAROSCOPIC/CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 31cm x 33cm Abdominal Fenestration
- Hook and Loop Tube and Line Control
- Fabric Reinforcement
- Armboard Covers

Code	Description	Packaging
89538	Laparoscopic/Cholecystectomy Drape, Sterile	Items/Case = 25
79538	Laparoscopic/Cholecystectomy Drape, Non-Sterile	Items/Case = 26

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15



ROBOTIC ABDOMINAL / PELVIC DRAPE, 264cm x 406cm

- 30cm x 10cm x 48cm Trapazoidal Abdominal Fenestration
- Clear Film Poly Pouches
- Tube Holders and Line Control
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
29283	Robotic Abdominal/Pelvic Drape, Sterile	Items/Case = 8
29284	Robotic Abdominal/Pelvic Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


MINOR PROCEDURE DRAPES



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Rectangular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Circular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Rectangular Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23



CIRCUMCISION DRAPE, 31cm x 51cm

- 2.5cm Circular Fenestration
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44972	Circumcision Drape, Sterile	Items/Case = 40

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS

THYROID DRAPES



THYROID DRAPE, 193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders


Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


UNIVERSAL PACKS

UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



Code	Description		Packaging
88759	Universal Pack – Major		Items/Case = 5
	Sterile Component Code	Non-Sterile Component Code	Contents
	89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 38cm x 56cm 1-Suture Bag
	34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
	34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
	34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm

UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES



Code	Description		Packaging
88769	Universal Pack – Minor		Items/Case = 14
	Sterile Component Code	Non-Sterile Component Code	Contents
	89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
		79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
		79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
		79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS



UNIVERSAL PACK I

Code	Description	Packaging
88878	Universal Pack – I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm
	79756	2-HALYARD BASICS* Side Drape, 193 x 112cm
	79754	1-HALYARD BASICS* Top Drape, 274 x 127cm
	79755	1-HALYARD BASICS* Bottom Drape, 193 x 173cm



UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
88774	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 2-Absorbent Towels, 38cm x 56cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS



U-BAR* PACK

Code	Description	Packaging	
88334	U-BAR* Pack	Items/Case = 16	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm	



U-BAR* PACK – IV

Code	Description	Packaging	
88773	U-BAR* Pack – IV	Items/Case = 10	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm	

HALYARD BASICS* U-BAR PACK I

Code	Description	Packaging	
88875	HALYARD BASICS* U-BAR Pack I	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89532	79532	1-HALYARD BASICS* U-Drape, 203cm x 315cm, with 10cm x 112cm Split	
89533	79533	1-HALYARD BASICS* Bar Drape, 244cm x 112cm	
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 112cm x 229cm	
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

GENERAL DRAPES & PACKS


UNIVERSAL DRAPES



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


GENERAL DRAPES & PACKS



TOP DRAPE, 276cm x 150cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70




SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

NEUROSURGICAL DRAPE




CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm

- CONTROL* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port

Code	Description	Packaging
89062	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS




C-SECTION DRAPES

C-SECTION DRAPE WITH CLEAR SCREEN, 264cm x 193cm x 305cm

- Triangular Fenestration
- 360° Fluid Collection Pouch
- Clear Screen
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44968	C-Section Drape with Clear Screen, Sterile	Items/Case = 6
44967	C-Section Drape with Clear Screen, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH POUCH AND FENESTRATION, 264cm x 193cm x 305cm

- 36cm x 31cm Adhesive Area with Fenestration
- 360° Fluid Collection Pouch
- CONTROL* Fabric Reinforcement
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44978	C-Section Drape with Pouch and Fenestration, Sterile	Items/Case = 6
44977	C-Section Drape with Pouch and Fenestration, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89528	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89529	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LAPAROSCOPY PACKS

LAPAROSCOPY PACK I



Code	Description	Packaging
88088	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89411		1-Under Buttocks Drape, Plastic, 75cm x 112cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89119	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm • 15cm x 20cm Rectangular Abdominal Fenestration • 7cm x 15cm Perineal Fenestration with Flap • Tube Holders



LAPAROSCOPY / PELVISCOPY PACK I

Code	Description	Packaging
88222	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm • 18cm x 23cm Trapezoidal Abdominal Fenestration • 8cm x 15cm Perineal Fenestration with Flap

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LAPAROSCOPY DRAPES



LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers


Code	Description	Packaging
89219	Laparoscopy/Pelviscopy Drape with Pouch, Sterile	Items/Case = 10



LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
89224	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10



LAPAROSCOPY "T" DRAPE, 249cm x 104cm x 251cm

- 15cm x 20cm Rectangular Abdominal Fenestration
- 7cm x 15cm Perineal Fenestration with Flap
- Tube Holders

Code	Description	Packaging
79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, Non-Sterile	Items/Case = 26

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY


DRAPES & PACKS

LITHOTOMY PACKS



LITHOTOMY PACK II

Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
		1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89101	79101	1-Half Drape, 152cm x 111cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LITHOTOMY DRAPES



LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings


Code	Description	Packaging
89203	Lithotomy Drape, Sterile	Items/Case = 16



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 80



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7cm Triangular Pouch

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


OPHTHALMIC DRAPE & PACKS



BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



OCUARTS* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders


Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS

ORTHOPAEDIC PACKS

ORTHO PACK I – MAJOR



Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
44682	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89331	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split 1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
89381	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



89345/79345

ORTHOARTS* HIP PACK WITH POCKETS

Code	Description	Packaging
88452	ORTHOARTS* Hip Pack with Pockets	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, in Overwrap, with Towel, Large
44682	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89821	79515	1-Stockinette, Large
89121	79121	1-Large Drape, 152cm x 193cm
89333	79333	1-Impervious Split Drape, Large, 152cm x 205cm with 10cm x 79cm Split
89345	79345	1-ORTHOARTS* Hip Drape with CONTROL* Plus Reinforcement, 289cm x 223cm x 320cm, with 21cm x 17cm Conformable Teardrop Fenestration 1-Yankauer Suction Tip 1-304cm Suction Tubing

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* LOWER EXTREMITY

Code	Description	Packaging
88471	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In 60cm x 120cm Overwrap, with Towel, Large
		1-Extremity Prep Pad
95221	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
44686	44668	1-AERO CHROME* Breathable High Performance Specialty Surgical Gown, X-Large, X-Long
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
		1-304cm Suction Tubing
89821	79515	1-Stockinette, Large
89331	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split
89278	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, with 6cm Conformable Circular Fenestration, 3-Drape Tape Strips
89053		3-Tape Strips, Sterile

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


EXTREMITY DRAPES



LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89016	Lower Extremity Drape, Sterile	Items/Case = 12



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



BILATERAL LIMB DRAPE, 193cm x 304cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders


Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBOARD COVERS, 295cm x 223cm x 366cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12



UPPER EXTREMITY DRAPE, 193cm x 305cm

- 3cm Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89271	Upper Extremity Drape, Sterile	Items/Case = 14
79273	Upper Extremity Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


HAND DRAPES



HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
77450	Hand Drape, Sterile	Items/Case = 11
13227	Hand Drape, Non-Sterile	Items/Case = 108



ORTHOARTS* HAND DRAPE, 193CM X 381CM


- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


HIP DRAPES



ORTHOARTS* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89345	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



ORTHOARTS* HIP DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89346	ORTHOARTS* Hip Drape, Sterile	Items/Case = 7
79346	ORTHOARTS* Hip Drape, Non-Sterile	Items/Case = 7



HIP DRAPE WITH CLEAR LEG POCKETS, 294cm x 223cm x 345cm

- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
89349	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS

KNEE ARTHROSCOPY DRAPES


KNEE ARTHROSCOPY DRAPE, 228cm x 315cm



- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 8cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89265	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150


KNEE ARTHROSCOPY DRAPE WITH ARMBOARD COVERS, 290cm x 223cm x 315cm



- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89371	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13


ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm



- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7

HALYARD BASICS* KNEE ARTHROSCOPY PACK



Code	Description	Packaging
88886	HALYARD BASICS* Knee Arthroscopy Pack	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88660	78660	1-HALYARD BASICS* Film Split Drape, 152cm x 178cm, with 10cm x 53cm Split
88661	78661	1-HALYARD BASICS* Knee Arthroscopy Drape with Fluid Collection Pouch, 229cm x 315cm
89564	79764	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 152cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10



HALYARD BASICS* KNEE ARTHROSCOPY DRAPE, 229cm x 315cm

- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
88661	Knee Arthroscopy Drape, Sterile	Items/Case = 20
78661	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric



Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders


Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric


Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12



HALYARD BASICS* SHOULDER SPLIT DRAPE, 285cm x 152cm

- 13cm x 34cm Split
- Fabric Reinforcement

Code	Description	Packaging
89559	Shoulder Split Drape, Sterile	Items/Case = 30
79759	Shoulder Split Drape, Non-Sterile	Items/Case = 50



HALYARD BASICS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 158cm

- Fluid Collection Pouch with Suction Port
- Tube Holders

Code	Description	Packaging
89558	Shoulder Arthroscopy Drape, Sterile	Items/Case = 30
79758	Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS

U & BAR DRAPES


U-DRAPE, 193cm x 305cm



- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190


LARGE U-DRAPE, 223cm x 317cm



- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15

ORTHOARTS* U-DRAPE, 223CM x 274CM



- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* BODY SPLIT DRAPE, 223CM x 294CM

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12



SPLIT DRAPE, 152cm x 178cm

- Impervious
- Polyethylene Material
- 10cm x 53cm Split


Code	Description	Packaging
89331	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



LARGE SPLIT DRAPE, 152cm x 205cm

- Impervious
- Polyethylene Material
- 11cm x 78cm Split

Code	Description	Packaging
89333	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100



HALYARD BASICS* BODY SPLIT DRAPE, 224cm x 295cm

- 10cm x 67cm Split
- Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89557	Body Split Drape, Sterile	Items/Case = 25
79757	Body Split Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


ORTHOPAEDIC DRAPES & PACKS



ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
36931	Ortho Split "T" Drape, Sterile	Items/Case = 10



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28




VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

ORTHOPAEDIC DRAPES & PACKS



ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
89093	Isolation Screen, Sterile	Items/Case = 24
79093	Isolation Screen, Non-Sterile	Items/Case = 30



STOCKINETTES

Code	Description	Packaging
89811	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Stockinette, Non-Sterile, Medium, 24cm x 117cm	Items/Case = 75
89821	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Stockinette, Non-Sterile, Large, 29cm x 117cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

UROLOGY DRAPES & PACKS




T.U.R. PACK

Code	Description	Packaging
77885	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component Code	Contents
		1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Absorbent Towel, 38cm x 50cm
76885		1-T.U.R. Drape, 289cm x 175cm x 170cm:
		• Attached Leggings
		• Fluid Collection Pouch with Screen, Suction Port and Tie
		• 54cm x 75cm Triangular Pouch
		• Finger Cot
		• 7cm Fenestration
		• 25cm x 20cm Fenestration



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
	Sterile Component Code	Non-Sterile Component Code
		Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm with 7cm x 15cm Perineal Fenestration and 23cm x 23cm Screen



PERCUTANEOUS NEPHROLITHOTOMY DRAPE (PCNL), 299cm x 182cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Fluid Collection Pouch
- 67cm x 61cm x 7cm Trapezoidal Pouch
- 18cm x 28cm Fenestration

Code	Description	Packaging
42427	Percutaneous Nephrolithotomy Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

STANDARD PACKS

ENT PACK

Code	Description	Packaging
ENTEPHYH-1A	ENT Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm 1-Bag, Ziplock, 16cm x 25cm 2-Blade, Surgical, Carbon Steel, 15 1-Bowl, 1000ml 1-Bowl, 250ml, Clear 1-Cover, Light Handle, Flexible, Green 2-Dish, Kidney, 700ml, Clear 1-Drape, Head Turban, 127cm x 91cm 1-Drape, EENT, Split, 193cm x 315cm 3-Gallipot, 150ml 5-Gauze, 8 Ply, 10cm x 10cm, XRD 1-Handle, Yankaeur, Bulb Tip, Non-Vented 1-Mayo Stand Cover, 58cm x 137cm 2-Needle, 18G x 38mm 2-Needle, 23G x 25mm 1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow 2-Syringe, 10ml, L/L 2-Syringe, 20ml, L/L 1-Tray, Sharps Transfer, Yellow 1-Tubing, Suction, Soft, 4m	

UNIVERSAL PACK

Code	Description	Packaging
GENUPHYH-1A	Universal Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm 1-Bag, Ziplock, 16cm x 25cm 1-Blade, Surgical, Carbon Steel, 15 1-Bowl, 1000ml 1-Cover, Light Handle, Flexible, Green 2-Dish, Kidney, 700ml, Clear 1-Drape, Bottom, 193cm x 195cm 2-Drape, Side, 195cm x 105cm 1-Drape, Top, 273 x 127cm 2-Gallipot, 150ml 5-Gauze, 8 Ply, 10cm x 10cm, XRD 5-Label, 125mm x 5cm 1-Marker, Skin, Standard Tip 1-Mayo Stand Cover, 58cm x 137cm 2-Needle, 23G x 25mm 1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow 1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm 1-Ruler, White, 15cm 5-Sponge, Abdominal, XRD, 30cm x 30cm 1-Syringe, 20ml, L/L 1-Tray, Sharps Transfer, Yellow	

NOTE: Please see package inserts for directions for use for each of these products..
 Always read the label and follow the directions for use.
 These products are not available for purchase by the general public.

STANDARD PACKS

CYSTO-GYNAE PACK

Code	Description	Packaging
OBGPHYH-1A	Cysto-Gynae Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bowl, 1000ml	
	1-Dish, Kidney, 700ml, Yellow	
	1-Dish, Kidney, 700ml, Clear	
	2-Drape, Leggings, 78cm x 122cm	
	1-Drape, Lithotomy T	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	1-Pad, Peri Winged 10cm x 27cm	

CAESAREAN PACK

Code	Description	Packaging
OBGCPHYH-1A	Caesarean Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bag, Suture, with tape, 20cm x 25cm	
	1-Bag, Zip Lock, 16cm x 25cm	
	1-Blade, Surgical, Carbon Steel, 23	
	1-Blanket, Baby, 75cm x 100cm	
	1-Bowl, 1000ml	
	1-Catheter, Y-Suction, 10FG x 45cm	
	1-Cover, Light Handle, Flexible, Green	
	2-Dish, Kidney, 700ml, Clear	
	1-Drape, C-Section, Fenestrated Incise	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	1-Handle, Yankauer, Bulb Tip, Non-Vented	
	1-Pad, Peri Winged 10cm x 27cm	
	1-Mayo Stand Cover, 58cm x 137cm	
	1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
	1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
	5-Sponge, Abdominal, XRD, 30cm x 30cm	
	1-Syringe, 10ml, L/S	
	1-Tray, Sharps Transfer, Yellow	
	1-Tubing, Suction, 3.6m	

NOTE: Please see package inserts for directions for use for each of these products..
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

STANDARD PACKS

LOWER EXTREMITY PACK

Code	Description	Packaging
ORTEPHYH-1A	Lower Extremity Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bag, Ziplock, 16cm x 25cm	
	1-Bandage, 15cm x 2.7m, Medium Weight	
	1-Bowl, 1000ml	
	1-Cover, Light Handle, Flexible, Green	
	2-Dish, Kidney, 700ml, Clear	
	1-Drape, Extremity, Extra Long, 289cm x 321cm	
	1-Drape, Split, 152cm x 178cm	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	5-Label, 125mm x 500mm, Blank	
	1-Marker, Skin, Standard Tip	
	1-Needle, 23G x 38mm	
	1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
	1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
	1-Ruler, White, 15cm	
	5-Sponge, Abdominal, XRD, 30cm x 30cm	
	1-Stockinette, Impervious, Large, 30cm x 122cm	
	1-Syringe, 20ml, L/L	
	1-Tray, Sharps Transfer, Yellow	

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

CUSTOM PROCEDURE PACKS (CPP)

THE FUTURE OF CUSTOM PROCEDURE PACKS FROM THE NAME YOU TRUST

For over 30 years HALYARD has been providing high quality surgical drapes and gowns, designed to deliver protection, comfort and optimal patient outcomes.



Leveraging Halyard's global expertise in supply chain, component sourcing, utilisation reviews and cost efficiencies, we are now extending our offering to provide ANZ clinicians with HALYARD branded Custom Procedure Packs (CPPs).

Based on years of working together with clinicians globally we know that some of the challenges you face in Operating Theatres include:

- Operational inefficiencies
- Lack of visibility to costs and inventories
- Suboptimal utilisation of Operating Theatre supplies
- Procedural supply waste
- Need for improved throughput
- Staff satisfaction and retention

That's where HALYARD comes in. When you make the decision to switch to HALYARD CPPs, our clinical experts work with you to ensure EVERY component is used EVERY time. It's part of our PACT pledge and the HALYARD CPP Solution.

NOTE: Please see package inserts for directions for use for each of these products..
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

INDEX

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
13227	40	42427	50	47619	15	79551	33	47622	16	89431	10
13231	46	42429	43	79117	31	79602	12	47630	16	89438	10
13233	38	42620	4	47624	16	79603	12	47628	16	89452	8
13294	46	42621	4	47631	16	79607	12	89203	33	89453	8
13296	7	43623	4	79147	16	79709	7	89204	33	89454	8
14007	12	43940	4	79204	33	79720	13	89206	21	89471	6
29219	38	43941	4	79206	21	79731	13	89207	21	89490	45
29220	38	44661	3	79207	21	79749	7	89208	21	89492	45
29240	19	44662	3	79208	21	79757	47	89209	21	89528	29
29241	19	44664	3	79209	21	79758	45	89219	31	89529	29
29258	12	44666	3	79224	31	79759	45	89224	31	89538	19
29259	12	44667	3	79233	19	79828	16	89229	19	89543	11
29269	12	44668	3	79240	20	79835	18	89233	19	89557	47
29283	20	44669	3	79255	22	79945	13	89236	17	89558	45
29284	20	44680	3	79265	42	88088	30	89241	20	89559	45
31445	18	44681	3	79267	42	88161	14	89255	22	89583	6
34901	26	44682	3	79273	39	88171	14	89265	42	89601	12
34902	26	44683	3	79275	47	88178	14	89267	42	89602	12
34903	26	44684	3	79279	39	88222	30	89271	39	89603	12
34904	26	44685	3	79282	20	88223	9	89272	40	89621	12
34905	26	44686	3	79304	35	88334	25	89274	45	89632	12
34906	26	44687	3	79330	47	88411	36	89275	47	89636	18
36930	43	44967	29	79333	47	88452	36	89276	38	89680	18
36931	48	44968	29	79335	17	88471	37	89278	38	89701	13
36932	44	44972	21	79343	16	88521	32	89279	39	89709	7
37498	12	44976	6	79345	41	88551	32	89282	20	89711	7
41725	4	44977	29	79346	41	88611	50	89291	39	89731	13
41727	4	44978	29	79349	41	88661	43	89301	46	89749	7
41739	4	70392	5	79351	10	88665	12	89304	35	89811	49
41740	4	72978	5	79356	35	88666	12	89311	46	89821	49
41741	4	72988	5	79360	48	88711	9	89331	47	89828	16
41870	4	76323	15	79363	46	88716	9	89333	47	89835	18
41871	4	76356	8	79368	6	88719	10	89343	16	89945	13
41872	4	76358	12	79371	42	88759	23	89345	41	90003	5
41873	4	77450	40	79408	33	88769	23	89346	41	90018	5
41874	4	77723	48	79411	34	88773	25	89347	46	90019	5
42217	12	77885	50	79414	34	88774	24	89349	41	90048	5
42218	12	78656	12	79415	34	88875	25	89352	10	90049	5
42224	12	78661	42	79430	10	88878	24	89353	15	90050	5
42225	12	78664	12	79452	8	88886	42	89356	35	90138	5
42231	12	79049	44	79453	8	89016	38	89368	6	ENTEPHYH-1A	51
42233	12	79061	44	79482	12	89023	15	89371	42	GENUPHYH-1A	51
42306	12	79070	44	79490	45	89053	13	89380	39	OBGPHYH-1A	52
42307	12	79075	28	79492	45	89062	28	89381	48	OBGCPHYH-1A	52
42309	12	79093	49	79505	49	89066	44	89408	33	ORTEPHYH-1A	53
42310	12	79095	27	79515	49	89070	44	89413	34	CPP	54
42366	12	79096	27	79528	29	89093	49	89414	34		
42367	12	79097	27	79529	29	47614	15	89415	34		
42426	44	47615	15	79538	19	89111	15	89421	11		



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



O&M Halyard Australia Pty Ltd.
Customer Care: 1800 664 227
or email CustomerService_ANZ@hyh.com

www.halyardhealth.com.au



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2022 O&M Halyard, Inc.
All Rights Reserved. 22-06:WFHY947v2